

# Prism C-380 Bariatric Ceiling Track Hoist



**PRISM MEDICAL UK**

*Moving, Handling and Bathing Solutions™*

# Prism C-380 Bariatric Ceiling Track Hoist

Bariatric clients often exceed a hoist's safe working load and this can create an unsafe working environment.

In cases where a client's weight, size or medical condition requires precise control when transferring or re-positioning, the Prism C-380 Bariatric Ceiling Track Hoist provides a safe and reliable solution.

The Prism C-380 is specifically designed to lift clients that exceed a standard hoist's maximum weight capacity.

As the strongest Ceiling Track Hoist in its class, the Prism C-380 safely transfers bariatric clients with respect and dignity.

**The Prism C-380 offers a safe working load of 380kg (just under 60 stone)**



## Key Features

- Safe working load of 380kg (just under 60 stone)
- Easy to operate pneumatic handset
- Emergency stop/emergency lower
- Manual or powered traverse available
- Smooth stop start to ensure user comfort
- Bariatric four point carry bar
- Carry bar fitted with Quick Release System for simple interchange of carry bars
- Easy to read LCD display
- Low profile hoist unit provides a great lift height

**To arrange a demonstration of the Prism C-380 please contact us on 0845 409 2000 or visit [www.prismmedical.co.uk](http://www.prismmedical.co.uk)**